



AcuPac Advanced[®]

BONE VOID FILLER ALLOGRAFT

AcuPac Advanced[®] is a **premier human tissue allograft consisting of frozen cancellous bone and demineralized cortical fibers** that function as a bone void filler to support bone repair.

BENEFITS	
Improved Alternative To Autograft	AcuPac Advanced [®] is a refined alternative to autograft, with the additional benefit of demineralized cortical fibers. It is a natural bone void filler surgeons can use to promote rapid, complete bone regeneration. AcuPac Advanced [®] also provides improved handling and osteoinductive payload.
Fill Gaps & Promote Bone Growth	The pre-milled, cancellous bone product combined with demineralized cortical fibers to fill the gaps in bone caused by surgical procedures, trauma, infection, or excision of tumor(s).
Frozen & Ready-To-Use	AcuPac Advanced [®] provides maximum surgical flexibility. After thawing, the bone product has a fibrous cancellous autograft consistency that can fill in defects of all shapes and sizes.
Rigorously Screened For Highest Quality Bone	AcuPac Advanced [®] lets you ensure the highest quality bone product for your patient. It goes through an exacting donor screening and testing process from an AATB-Accredited Tissue Bank.

INDICATIONS FOR USE

AcuPac Advanced® is an **allogeneic bone graft, which is used as bone void filler for bridging gaps in bone caused by surgical procedures, trauma, infection or excision of tumor(s)**. The allograft can be used in a number of orthopedic, spine, general and reconstructive surgical applications. AcuPac Advanced® is **intended for one-time use only, for a single patient**, and by a licensed physician, dentist or podiatrist.

AcuPac Advanced® combines pre-ground cancellous of 1mm – 4mm in size, with the addition of demineralized cortical fibers that is aseptically processed. **Because of minimal processing, AcuPac Advanced® retains more of the natural components found in bone.**

AcuPac Advanced® is **supplied frozen and packaged sterile in a peel pouch contained within another peel pouch**. Allograft volume is indicated on the package label.



AcuPac Advanced®

BONE VOID FILLER ALLOGRAFT

90-F1040002

2.5 cc

90-F1040010

10 cc

90-F1040025

25 cc