



AcuPure™ DBM Putty

BONE VOID FILLER STERILE ALLOGRAFT

AcuPure™ DBM Putty is a **gamma-irradiated, sterile human tissue allograft consisting of 100% human demineralized bone allograft**. These allografts function as a bone void filler to support bone repair.

BENEFITS

Promotes Bone Formation

Innovation often follows the straightest path. In bone grafting, attention to basic fundamentals has always been the key to the best and most predictable clinical outcomes. Historically, bone itself is the best substrate for new bone formation. This is the philosophy behind AcuPure™ DBM Putty, which is comprised of 100% human bone.

Flowable, Formable Putty

AcuPure™ DBM Putty shapes and conforms to fill unique bony voids and has excellent handling characteristics to retain its shape without migrating or washing away with lavage.

Sterile & Ready-To-Use

AcuPure™ DBM Putty is gamma irradiated and stored at ambient room temperature. No prep time or reconstitution required as it is ready for immediate use.

Rigorously Screened For Highest Quality Bone

AcuPure™ DBM Putty lets you ensure the highest quality bone product for your patient. It goes through an exacting donor screening and testing process from an AATB-Accredited Tissue Bank.

INDICATIONS FOR USE

AcuPure™ DBM Putty is an allogeneic bone graft, which is used as bone void filler for bridging gaps in bone caused by surgical procedures, trauma, infection or excision of tumor(s). The allograft can be used in a number of orthopedic, spine, general and reconstructive surgical applications. AcuPure™ DBM Putty is **intended for one-time use only, for a single patient**, and by a licensed physician, dentist or podiatrist.

AcuPure™ DBM Putty is **supplied at ambient temperature and packaged sterile in a peel pouch contained within another peel pouch**. Allograft volume is indicated on the package label.



AcuPure™ DBM Putty

BONE VOID FILLER STERILE ALLOGRAFT

PRODUCT CODE	VOLUME
90-00003301	1 cc
90-00003302	2.5 cc
90-00003305	5 cc
90-00003310	10 cc